

C E R T I F I C A T E O F C O M P L I A N C E

LOT #: 602502

KEY SCIENTIFIC PRODUCTS  
 1113 E REYNOLDS STREET  
 STAMFORD, TEXAS 79553  
 1-800-THE-1-KEY

THIS LOT HAS BEEN TESTED ACCORDING TO THE CRITERIA LISTED BELOW AND MEETS THE PERFORMANCE STANDARDS AS ESTABLISHED BY KEY SCIENTIFIC PRODUCTS TO MEET NCCLS GUIDELINES.

PRODUCT: BUTYRATE DISCS  
 CATALOG #: K9131B

LOT #: 602502

PACKAGED DATE = 2021/07/30 EXP DATE = 2023/07/30

CONTROL LOT# 594900

TESTED BY: BLC

		CURRENT POSITIVE	CURRENT NEGATIVE	CONTROL POSITIVE	CONTROL NEGATIVE
25923	STAPHYLOCOCCUS AUREU	BLUE		BLUE	
49143	MORAXELLA CATARRHALI	BLUE		BLUE	
43069	NEISSERIA GONORRHOEA		COLORLESS		COLORLESS
13813	STREPTOCOCCUS AGALAC		COLORLESS		COLORLESS
6305	STREPTOCOCCUS PNEUMO		COLORLESS		COLORLESS
19615	STREPTOCOCCUS PYOGEN		COLORLESS		COLORLESS
27823	STREP CONSTELLATUS		COLORLESS		COLORLESS

DRY PAPER COLOR - WHITE  
 POSITIVE IS BLUE AFTER INCUBATION OF 15 MINUTES  
 NEGATIVE REMAINS COLORLESS