

C E R T I F I C A T E   O F   C O M P L I A N C E

LOT #: 639911

KEY SCIENTIFIC PRODUCTS  
1113 E REYNOLDS STREET  
STAMFORD, TEXAS 79553  
1-800-THE-1-KEY

THIS LOT HAS BEEN TESTED ACCORDING TO THE CRITERIA LISTED  
BELOW AND MEETS THE PERFORMANCE STANDARDS AS ESTABLISHED  
BY KEY SCIENTIFIC PRODUCTS TO MEET NCCLS GUIDELINES.

PRODUCT: HIPPURATE DISCS 15/BTL  
CATALOG #: K710315

LOT #: 639911

PACKAGED DATE = 2026/06/17   EXP DATE = 2028/06/17

CONTROL LOT# 580800

TESTED BY: BLC

|       |                      | POSITIVE<br>CURRENT | NEGATIVE<br>CURRENT | POSITIVE<br>CONTROL | NEGATIVE<br>CONTROL |
|-------|----------------------|---------------------|---------------------|---------------------|---------------------|
| 12386 | GROUP B STREP        | +                   |                     | +                   |                     |
| 19615 | STREPTOCOCCUS PYOGEN |                     | -                   |                     | -                   |

REAGENT # 614100  
DISK COLOR WHITE   INCUBATE FOR 2 HOURS  
POSITIVE: PURPLE COLOR AFTER ADDITION OF DEVELOPER.  
NEGATIVE: NO COLOR CHANGE AFTER ADDITION OF DEVELOPER.  
CHECK HIPPURATE VOL /DISK  
WEIGH 10 BLANK DISKS OF SAME SIZE  
WEIGH 10 FINISHED HIPPURATE DISKS  
DIFFERENCE SHD BE 1.8 TO 2.2 MG PER DISK